



Automated Mixing Device (MolGeniusComposer)



Modular pipette (5mL)



Component layout diagram



Multi-channel magnetic stirrer with individual speed control

	tube 0	tube 1	tube 2	tube 3	...
round 0	L:3.0	L:3.0	L:3.0	L:3.0	...
round 1	d:1.0	d:1.0	a:1.0	a:1.0	...
round 2	c:1.0	e:1.0	e:1.0	e:1.0	...
round 3	a:1.0	a:1.0	b:1.0	d:1.0	...
round 4	b:1.0	c:1.0	d:1.0	c:1.0	...
round 5	e:1.0	b:1.0	e:1.0	b:1.0	...


A55 SolutionMixingProtocol.csv (Partial View)

- Read tasks from the CSV protocol file and add them to the task pool.

Priority	From	To	Volume
0	Bottle L	Tube 0	3 mL
1	Bottle d	Tube 0	1 mL
2	Bottle c	Tube 0	1 mL
3	Bottle a	Tube 0	1 mL
4	Bottle b	Tube 0	1 mL

Task Pool

- Sort tasks by priority and merge them according to the "from" field.



- Execute the addition tasks.